

TiFast Titanium

• Dr. Ivako Rossi.

• billetter: 135mm

• TiFast VAK Rebody laminar Sinter.

• bars → 7 to 60 mm finished bars, not drawn wires

• wires → 7.8 (Noll diameter) → 13 mm

bars 13 mm to 62 mm

peeling in bars 8 to 60 mm

bill poling, US, Die Restant. → 8 inches

peeling 16 mm on Turnage → 125 mm on cable
laminar.

ISO 9001 - ISO 9100, WADCAP System and Heat Treatment.

US 2 level 3, but per Wadcap Environments.

Challenge:

TiFast Titanium

• ⇒ Collocatio

• ELG →

• PCC → market preparation

• Good Runnings: Tings

• Small Runnings: billets + peeling + Oil ⇒ sell to

processor

• 0.16 / O₂:

Tings ELI 0.16 → 0.11

base of selling: leg

entry:

billets peeling

blanks coils - interest for long term markets,
small wires

medical Tings 7 to 60 ~~mm~~

Fasteners